

584

THE STUDY OF NMR IMAGING IN ANIMAL MODEL.
T.Yamane,H.Ikehira,N.Fukuda,S.Torii,
S.Shibata and Y.Tateno. Natinal Institute
of Radiological Sciences,Chiba

These animal models were produced for dynamic function study of kidney by NMR renogram with Gd-DTPA which is contrast enhancement agent for NMR imaging.

The animals used that left ureter occlusion, right renal artery stenosis and hazard by cisplatin of rabbit.

The determination of spin relaxation time (T1) were measured using Asahi Mark-J(0.1T resonance frequency for hydrogen of 4.5MHz, the receiver coil was 24cm in diameter.

These results, in the normal kidney Gd-DTPA was 50% decay about 20 min. and 90% decay about 90 min. after administration of dose. On the other hand in the disease kidney decay curves were shows different then normal kidney. that in the ureter occlusion model shows correction of Gd-DTPA is long time and in the renal artery stenosis model shows excretion was very slow to kidney from blood vessel, and in the cisplatin administrated rabbit was not excrete of Gd-DTPA.

We confirmed these examination animal disease model are useful for dynamic function study by NMR-renogram with Gd-DTPA.

585

THE EFFECT OF CONTRAST ENHANCEMENT IN RABBIT TUMOR BY GD-DTPA ON MRI.S.Nawano.
National Sakura Hospital.Chiba.N.Arimizu,
T.Miyoshi,M.Saito,M.Ozaki.Chiba University
School of Medicine,Chiba.

Magnetic resonance imaging was performed on rabbits with VX-2 tumors, before and after intravenous administration of gadolinium-DTPA(Gd-DTPA) MR images were obtained with a VISTA-MR (Picker),superconducting magnet system operating at 0.26 Tesla.Contrast-enhanced CT of tumors were performed to compare with MRI.The IR(1500/500/40) image and SE (300/24) image were obtained,before and after Gd-DTPA.The capsule of VX-2 tumor intensity was increase not only on MRI but on X-CT.The intensity of necrosis was increase very slowly on MRI. After Gd-DTPA IR images had good contrast between necrosis and not. We thought Gd-DTPA was good contrast agent for MRI.

587

NMR:T1,T2-SENSITIVITY MAP OF SIGNAL INTENSITY TO DIAGNOSIS OF NMR-CT
S.Torii,H.Ikehira,N.Fukuda,Y.Tateno,
M.Endou,T.Matsumoto,T.Iinuma, Division of
Clinical Research National Institute of
Radiological Sciences,Chiba, Y.Ueshima,
J.Nishizawa, Asahi Medical Co.,Ltd.

In this study, we prepared T1,T2-sensitivity (total relaxation sensitivity) map as a means of presentation of the correlation between T1,T2 change and signal intensity change of various tissues. Ratio between fractional rate changes of each relaxation time against entire relaxation time and those of signal intensity against entire signal intensity, are measured respectively by means of logarithmic differentiation, by which the sensitivity of each T1 and T2 was calculated. The summation of respective relaxation sensitivity was designated as total relaxation sensitivity;(S.total). The resulting S.total is classified into seven different color indices according to its value. The classified S.total are plotted on T1 axis and T2 axis. The classified S.total express generally the degree of dependence of respective relaxation time and we discuss about the precaution and limitation of T1,T2-sensitivity map application.

588

IN VIVO MRI AND P-31 NMR SPECTROSCOPY IN EXPERIMENTAL HEAD INJURY.
T. Hashimoto, L.H. Pitts. Neurosurgery
Chiba Cancer Center Hospital, Chiba.
Neurosurgery, UCSF, USA.

The pathophysiological consequences immediately after head injury are complicated. The accurate assessment of brain damage is indispensable in the management of head injury victims, though methodologically difficult in clinical settings.

We used in vivo MRI (resistive, 0.35T) and P-31 NMR spectroscopy (5.6T) to evaluate the development of brain edema and the changes in brain phosphorus metabolism following mechanical head trauma in the rat. MRI 2 hours after injury readily demonstrated brain edema as high signal intensity region at the impacted cortex (spin echo pulse sequence, TR 2sec, TE 28msec). P-31 NMR spectroscopy showed the decrease of the PCr/Pi ratio and lower intracellular pH according to the magnitude of the impact pressure.

We conclude that NMR techniques will be of quite use in the diagnosis and management of brain edema and cerebral metabolic deterioration associated with mechanical head injury.